

Disposable Head Blocks - Foam Orange Blocks

Morrison Medical's Head Blocks are a foam-filled, light-weight head immobilizer. This device is affordably priced to be disposable, yet is durable enough to be reused. The highly visible orange vinyl exterior is both puncture and water resistant, while offering a writable surface for easy ownership and patient identification. Morrison Medical offers compatible disposable head and chin straps that are ideal for use with Head Blocks.

- Rectangular (4 x 5 x 10") foam blocks sealed in highly visible orange vinyl, complete with handle and writable surface
- Offers firm yet lightweight support.
- X-ray translucent.
- Tape can be used to secure patient, but works great with Morrison Medical Disposable Head/Chin Straps (Item 1475).
- Adjustable to fit children and adults.
- Compatible with helmets.
- Easy to use and quick application.
- Disposable.
- Puncture and water resistant.
- Head Blocks float - ideal for water rescue.



Disposable Head Blocks

Item 0320:

Head Blocks, pair

Item 0320-9:

Head Blocks, case of 9 pair

Item 0320-30:

Head Blocks, case of 30 pair

Disposable Head Blocks with Straps

Item 0322:

Head Blocks, pair (0320) and disposable head and chin straps (1475, 2)

Item 0322-9:

Case of 9 sets of item 0322

Item 0322-30:


Case of 30 sets of item 0322

Disposable Head/Chin Strap Self Adhesive

Item 1475, Each

Item 1475-100, Case of 100

- Disposable self adhesive head/chin strap, 2x36" .
- Easy and quick application - peel backing off, set in place and adhere to self to lock in place.
- Cut to remove.
- Available in blue and orange.



*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*