

Sticky Blocks Head Immobilizer

Morrison Medical's Sticky Block Head Immobilizer is constructed of foam contoured for maximum head stabilization and patient comfort. Adhesive backing of blocks allows for a quick and firm application and leaves no residue on backboards after removal. Adjustable for patients of all sizes.


Sticky Blocks Head Immobilizer

Please Call for Pricing

*Item 1400:
Sticky Blocks, pair*

*Item 1400-12:
Sticky Blocks, case of 12 pair*

*Item 1400-50:
Sticky Blocks, case of 50 pair*



- Contoured for maximum head stabilization
- X-ray translucent
- Universal size (3" x 5 1/2" x 10 3/8")
- Adjustable to fit children and adults; compatible with helmets
- Comfortable medium density foam
- Provides firm support
- Easy to use and quick application with adhesive backing
- No sticky residue after removal
- Disposable
- Works great with Morrison Medical Disposable Head & Chin Straps (Item 1475)

*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*

