


Morrison Medical's Big Blue Head Immobilizer

- Vinyl coated foam for easy disinfecting.
- Extra large pillows provide maximum support.
- Easy access to ears for patient assessment.
- Pillows adjust and firmly attach to base with wide hook and loop.
- Secures to backboard with mounted straps and D-rings.
- Includes base with attached straps, two large pillows and two head/chin straps.
- Adaptable for pediatric patients with the use of pediatric head blocks (#1432).
- Available in blue, orange and black.




Big Blue Head Immobilizer

Please Call for Pricing

1420	Big Blue Head Immobilizer, Each
1420-6	Big Blue Head Immobilizer, Case of 6
1422	Big Blue Pillows, Pair
1422-6	Big Blue Pillows, 6 Pair
1423	Big Blue Head & Chin Straps, Pair
1424	Big Blue Base Pad, Each
1432	Blue Pediatric Pillows, Pair

Pediatric Head Immobilizer Blocks

- Compatible with Big Blue and most other popular head immobilizers.
- Vinyl coated foam for easy disinfecting.
- Easy access to ears for patient assessment.
- Pillows adjust and firmly attach to base with hook and loop.
- 2½" x 5½" x 3¾" pillows.
- Available in blue.



*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*