

The Perry Pouch I & II From Morrison Medical

The Perry Pouch *Please Call for Pricing*

Item #3000:
Perry Pouch I, single handles, each

Item #3001:
Perry Pouch II, double handles, each



The Perry Pouch II features double handles for optimal patient handling.



The Perry Pouch I & II cinch around the patient for a custom fit and maximum stability during transport.

The Morrison Medical Perry Pouch was designed for fast and safe patient transport. When a patient is lifted, The Perry Pouch cinches around them for a custom fit and maximum stability. Ideal for transporting patients up to 500 pounds through confined areas and stairways, The Perry Pouch is compact and lightweight at 1 pound, 4 ounces and comes with a vinyl carrying case for easy storage. The orange vinyl surface and 1" black polypropylene webbing provide a durable surface that cleans quickly and easily. Every Perry Pouch comes with no-slip grip, rubberized and pre-molded handles to provide comfort and prevent fatigue to those transporting the patient. The Perry Pouch II features double handles at each of the four carrying points for even more "hands-off" control.

*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*

