

Disposable Polypropylene Straps

Morrison Medical disposable straps secure a patient to prevent further injury during transport. Economically priced and designed to be disposable, but durable enough to clean and reuse. Straps are easily adjusted for a custom fit and constructed of 2" polypropylene webbing with an 1800 pound break strength. Morrison Medical disposable straps meet or exceed FDA requirements and KKK Specification 3.11.8.1.



Disposable Poly Strap, Item 1392



- Available in white (WH), and orange (OR).
- 1,800 pound break strength

Disposable Poly Straps with Plastic Side Release Buckle

- 1390** 2 piece, 5', Loop Lok Ends, each
- 1390-SET** 2 piece, 5', Loop Lok Ends, set of 3
- 1390-36** 2 piece, 5', Loop Lok Ends, case of 36
- 1392** 1 piece, 9', each
- 1392-SET** 1 piece, 9', set of 3
- 1392-36** 1 piece, 9', case of 36
- 1394** 2 piece, 5', Plastic Swivel Speed Clip Ends, each
- 1394-SET** 2 piece, 5', Plastic Swivel Speed Clip Ends, set of 3
- 1394-36** 2 piece, 5', Plastic Swivel Speed Clip Ends, case of 36



Disposable Poly Strap, Item 1394

Disposable Poly Straps with Double Adjust Side Release Buckle

- 1396** 1 Piece, 6', each
- 1397** Replacement strap for item 1396 ONLY, 6', each
- 1396L** 2 Piece, Loop-lok ends, 7', each
- 1397L** Replacement strap for item 1396L ONLY, 7', each

Additional strap combinations and custom straps available. Call 800.438.6677 for more information.

Screen printing available. Make your straps easily identifiable with screen printed contact information or your logo. Call 800.438.6677 for more information.



Disposable Poly Strap, Item 1390

*For more information call 800-438-6677,
or visit www.MorrisonMed.com.*