


Foam EEZ Head Stabilization System

The Foam EEZ is an economical, all weather, disposable head immobilization system constructed of comfortable, yet sturdy, foam. The specially designed foam pillows feature easy access to ears for patient assessment. The system includes two foam pillows, a pair of black hook and loop head and chin straps and four cervical alignment foam pads. Foam EEZ Head Blocks are also available in pairs separate from the full system and work great with Morrison Medical Disposable Head & Chin Straps (item 1475). Compatible with all backboards. X-ray translucent. Latex-free.

Dimensions: Foam EEZ Head Blocks, 8.75" x 4" x 4.75"; cervical alignment foam pads, 4½" x 4½"; head and chin straps, 41.75" extended, loop section 24.5", hook sections 8.75" (at each end).

- Adjustable to fit children and adults
- Compatible with helmets
- Comfortable medium density foam
- Provides firm and comfortable support
- Easy to use and quick application on any backboard
- Disposable
- X-ray translucent
- Latex-free




Foam EEZ Head Stabilization System, Item 1401

Foam EEZ

Please Call for Pricing

ITEM	DESCRIPTION
1401	Foam EEZ Head Stabilization System, Each
1401-10	Foam EEZ Head Stabilization System, Case of 10 Sets
1402	Foam EEZ Head Blocks, Pair
1402-30	Foam EEZ Head Blocks, Case of 30 Pair



Foam EEZ Head Blocks, Item 1402

A Great Accessory for the Foam EEZ Head Blocks

Disposable Head/Chin Strap Self Adhesive

Item 1475, each



- Disposable self adhesive head/chin strap, 2 x 36".
- Easy and quick application - peel backing off, set in place and adhere to self to lock in place.
- Cut to remove.
- Available in blue and orange.

*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*