



3735 Paragon Drive
Columbus, Ohio 43228
800-438-6677 * 877-469-9696 Fax
www.morrisonmed.com

Perry Pouch Instructions for use

Description/Intended Use:

The Morrison Medical Perry Pouch is a soft stretcher that cinches around patient when lifted. It is designed for fast, safe patient transport.

Instructions:

Place Perry Pouch under patient. Lift patient with handles.
Supports up to 500 pounds.

Only qualified and properly trained personnel should use the Morrison Medical Perry Pouch. Any persons using the Perry Pouch should practice application and be thoroughly familiar with all medical protocols and aspects of its use **BEFORE** using. These instructions are provided as a guide only. Use protocols and techniques established by your local medical control.