

The Quick Transfer (QT)

The QT is easy to use, affordable, completely disposable and an effective way to transport patients who do not require complete spinal immobilization. The product is constructed of 4 mil blue plastic corrugated material. The QT comes complete with a built in head immobilization unit and includes 3 disposable stick & cut straps and 1 disposable head strap.



What makes The QT the superior choice?

- Folds down to 21" x 7" for easy storage
- Converts to pediatric size
- Tapered at the bottom
- Durable when wet



Item # 3045

Also available in case of 10.

*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*