



3735 PARAGON DRIVE  
COLUMBUS, OHIO 43228  
800.438.6677 \* 614.469.9696 FAX  
WWW.MORRISONMED.COM

## **Soft Stretcher Instructions for Use**

### **Description / Intended Use**

Morrison Medicals Soft Stretchers are specially designed for moving bariatric patients. It is great for transporting patients out of confined areas or stairways. Each size supports up to 2000 pounds.

### **Instructions**

Place Soft Stretcher under patient. Lift patient gently with handles. Not intended for use with suspected spinal injured patients.

Only qualified and properly trained personnel should use Morrison Medicals Soft Stretchers. Any persons using the Soft Stretchers should practice application and be thoroughly familiar with all medical protocols and aspects of its use BEFORE using. These instructions are provided as a guide only. Use protocols and techniques established by your local medical control.